

## ADME: Cytochrome P450 inhibition in human liver microsomes

### Background:

Co-administration of drugs can result in drug-drug interaction, as drugs compete for the same enzymes. Inhibition of CYP450 enzymes is a principal mechanism of metabolism-based drug-drug interactions and a common cause of adverse drug events. Therefore, one of the crucial properties in early drug discovery is the assessment of the potential of test compounds to inhibit a specific CYP450 isoform. This data is useful for further investigations of clinical drug-drug interactions (DDI).

The human liver microsomal assay is accepted as the 'gold standard' for *in vitro* DDI assessment as it is closest to the native enzyme environment. This assay is therefore used to characterise development compounds where data will be submitted to regulatory authorities.

### Assay description

#### CYP450 isoforms

CYP1A2, CYP2C9, CYP2C19, CYP2D6, CYP3A4

#### Compound concentration

10µM (single point) or 0-100µM (IC<sub>50</sub>)

#### Compound requirements

1-2 mg of dry matter

#### Incubation details

isoform specific substrate (see Table 1)

isoform specific time of incubation

#### Assay controls

known isoform specific positive control

(see Table 1 and Figure 1)


#### Detection method

LC-MS/MS with internal standard

**Table 1.** CYP450 isoform specific substrates and positive controls

Isoform	Substrate	Positive control
CYP1A2	Phenacetin	Fluvoxamine
CYP2C9	Diclofenac	Sulfaphenazole
CYP2C19	S-mephenytoin	Ticlopidine
CYP2D6	Bufuralol	Quinidine
CYP3A4	Midazolam*	Ketoconazole

\* Testosterone also used as substrate



**Figure 1.** Comparison of IC<sub>50</sub> values for CYP450 isoform specific positive control with literature values<sup>1,2,3</sup>

<sup>1</sup> Obach 2006, J Pharmacol Exp Ther 316, 336

<sup>2</sup> Zientek et al. 2008, J Pharmacol Toxicol Method 58, 206

<sup>3</sup> Zamboni et al. 2010, Drug Metabol Lett 4, 120

**Assay details adjustable to client's and/or project specific requests**

Contact us:

Adrijana Vinter, Business Development  
adrijana.vinter@glpg.com  
+385 91 265 5527

Mila Vrančić, Business Development  
mila.vrancic@glpg.com  
+385 91 265 5528