

## ADME: Permeability- Caco-2

### Background:

Caco-2 cell line is a line of human epithelial colon adeno-carcinoma cells. When cultured under specific conditions they become differentiated and polarised such that morphologically and functionally resemble the epithelial cells of human small intestine. Therefore, they are considered to be a golden standard for *in vitro* prediction of permeability and for estimation of the potential for different efflux/uptake transporters (most importantly P-glycoprotein, P-gp) to transport the compound of interest.

Transport of compound through Caco-2 cell monolayer is used for measuring the rate of membrane transport, expressed as apparent ( $P_{app}$ ) permeability. Resulting  $P_{app}$  values from both transport directions (apical to basolateral side and basolateral to apical side) are used for calculation of efflux ratio. In addition, the inclusion of the P-gp inhibitor could identify whether membrane transport is mediated by P-glycoprotein.

### Assay description

#### Cells

Caco-2

#### Direction

apical to basolateral (A2B)

basolateral to apical (B2A)

with or without P-gp inhibitor (elacridar, verapamil)

#### Compound concentration

10 $\mu$ M (1% DMSO)

#### Compound requirements

50 $\mu$ l of 10mM stock solution or

1-2 mg of dry matter

#### Incubation details

medium: Dulbecco's PBS (pH 7.4)

calibration curve: optional

#### Detection method

LC-MS/MS with internal standard

#### Results

$P_{app}$ (A2B),  $P_{app}$ (B2A)

efflux ratio, recovery

Caco-2 cells seeding on  
24-well plates  
(21 days prior experiment)



incubation of test compound at  
37°C for 1 hour  
(n=2)



sampling at 0min and 1h



quantification LC-MS/MS

## **ADME: Permeability- Caco-2**

### **Assay controls**

reference compounds: amprenavir and propranolol  
cell monolayer integrity control: lucifer yellow

**Assay details adjustable to client's and/or project specific requests**